DIFFUSION STARTER KIT

the standard for your lab.

PermeGear is pleased to announce our

Diffusion Starter Kit

This is intended for anyone new to diffusion testing and is all you require to begin. We are offering this package at a discounted rate for students and those who are just starting out.



or for 2 Standard Franz Cell.



This Side bi Side Cell has the common 9mm orifice diameter with the 3.4ml volumes.

The Starter Kit includes:

- 1 ea. **H1C Stirrer with Clamp**
- 1 ea. 9mm Side bi Side Cell with 3.4ml volume, stirbars and 7/25 stoppers
- 2 ea. 9mm jacketed Franz Cell with flat ground joint, 5ml receptor volume clear glass, stirbar and clamp
- 2 ea. Cell Stand for jacketed Franz Cell
- 1 ea. Sample Pack of Pippette Tip

For more informations please email us to: support@permegear.de



This diffusion cell is the most common Franz Cells. It is a jacketed Franz Cell with 9mm orifice diamter, flat ground joint, and 5ml receptor volume. It is the defacto standard for working with human skin and may be used with any membrane4mm or less thickness.

INCLUDED

Assemblies consisting of a three piece flexible assembly, the first designed specifically for insertion into the center of Franz Cells for more accurate sampling.

SES GmbH - Analytical Systems

Friedhofstrasse 7-9 - 55234 Bechenheim / Germany Tel: ++49(0)6736 1301 - E-Mail: ses_analysesysteme@t-online.de Visit our Webside: www.permegear.de or www.franzzellen.de

